

Model: CMRTSG400

400L Glass Door RTS Cabinet

Overview

- +15°C to +25°C temperature range
- Core temperature display, represents your vaccine vial through **iCOOL** temperature probe
- External digital controller and temperature display
- Battery Back up
- Annual service & Calibration available
- Fitted lock with 2 keys

5 YEARS PARTS & LABOUR WARRANTY

ANNUAL SERVICE & CALIBRATION AVAILABLE



Key Features



Controller

Controller system. Read, record, reset with the touch of one button. Adjustable temperature setting.



Battery

Backup should there ever be a power outage the battery backup will continue to monitor the internal temperature of the freezer.



Porthole

Factory fitted porthole enables quick and easy access for any secondary temperature monitoring probes.



LED Lighting

Low voltage LED lighting improves running cost as no heat is produced. Long lasting.



Flexible storage

Multiple shelf locations. Back stop on every shelf provides a flexible storage solution.

Details

Model: CMRTSG400

| | |
|--|--|
| Gross Capacity (litres) | 408 |
| Net Capacity (litres) | 400 |
| Unit External Width (mm) | 595 |
| Unit External Depth (mm) | 695 |
| Unit External Height (mm) | 1845 |
| Product Weight (kg) | 84 |
| Packed Width (mm) | 650 |
| Packed Depth (mm) | 750 |
| Packed Height (mm) | 1910 |
| Packed Weight (Kg) | 92 |
| Temperature Range | +15°C to +25°C |
| External Finish, Internal Finish | Painted Steel, Painted Steel |
| Battery Back Up | Yes |
| Door Finish, Door Type, Door Style | Painted Steel, Glass, Hinged, Right Hand |
| Open Door Sensor | Yes |
| Lock | Yes |
| Adjustable Feet (number) | Yes, 2 |
| Porthole | Yes |
| Porthole Location | Rear, Middle |
| Temperature Display on Panel, Temp Alarm, Remote Alarm Contacts, Audible Alarms, Visual Alarms | Yes |

| | |
|-----------------------------------|-------------------------------------|
| Number of Shelves | 7 (6 full size, 1 small base shelf) |
| Number of Shelf Positions | 10+ |
| Shelf Material | Plastic Coated Wire |
| Interior Light | Yes |
| Light Type | LED |
| Light Activation | Switch |
| Fan | Yes |
| Spark Free | No |
| Controller Type | Digital, External |
| Number of Probes | 2 |
| Probe Type | Air & Load |
| Ambient Temperature Range | +10°C to +30°C |
| Refrigerment Type | R600a |
| CFC Free, HCFC Free, Ammonia Free | Yes |
| Voltage Supply | 220-240V/ 50Hz |
| Certification | CE/UK |
| Noise Level (dbA) | 50 |
| Input Power (W) | 240 |
| Cable Length (Cm) | 20 |
| Energy Consumption (KWh/annum) | 2.6 |

Related Products

**CMRTSS400**400L Solid Door
RTS Cabinet**CMDL01**Advanced Data Logger
with LCD Display**COMBI01**Keyless
Combination Lock